



GE BRIGHTSPEED ELITE 16 CT1

System ID:	CT1
Manufacturer:	General Electric
Model:	BrightSpeed Elite 16
Year of production:	2011 / Tube YOM - 2018
No. of slices:	16
Software:	13BW29.5
Gantry opening:	60cm
Max table load:	205kg
Options:	Smart Prep, Smart Speed, Helical Tilt, AutoMA, Power440, Patient-16-slice 3000 Image Series, Direct-3D, 90kV, A, CopyComposer, Thin Twin Helical, Direct-MPR, Data Export, Connect Pro, Volume Viewer, EnhancedRecon, ASIR, Real Time CT Fluoro, Enhanced Xstream Injector

ELECTRICITY

Power cable length
15m - 20m

Maximum power parameters
120kW

Type of socket required for power supply
fixed socket 5P 125A 400V IP67

Cable length up to 50m
**5 x 35mm² (copper cable)
or 5 x 65mm² (aluminum cable)**

External connection protection
160A

TRAILER DETAILS

Type	Length
CALUMET	13.58m
Min. area for placement	Width
15.50m x 6.50m	2.46m
	4.78m (unfolded)
Weight	Height
30 200kg	3.95m

UNIT FEATURES

- Air-conditioning and recuperation system in all rooms
- Storage cabinets
- Sound system and FM tuner
- Hydraulic patient lift
- Sink with running water
- IT/Network connections in operator room
- Patient changing room

ADDITIONAL FEATURES

- Reference clinician workstation
- CD/DVD burner
- RIS/PACS/Network connection
- Injector: Medrad Stellant - dualhead

MAINTENANCE AND SERVICE

Maintenance and service contracts included