



## GE OPTIMA CT 520 CT4

<b>System ID:</b>	CT4
<b>Manufacturer:</b>	General Electric
<b>Model:</b>	Optima 520
<b>Year of production:</b>	2013 / Tube YOM - 2020
<b>No. of slices:</b>	16
<b>Software:</b>	12BW31.11_SPI-3-1
<b>Gantry opening:</b>	70cm
<b>Gantry serial number:</b>	354316HM2
<b>Max table load:</b>	227kg
<b>Options:</b>	Smart Prep, Helical Tilt, AutomA, Patient-16-slice 3000 image series, Direct-3D CopyCompressor, Direct-MPR Data Export, Lite-Power 350, Helical-120-sec, 0.8 Second Scan, VolumeViewer, Connect Pro, EnhancedRecon, ASIR

### ELECTRICITY

Power cable length  
**15m - 20m**

Maximum power parameters  
**120kW**

Type of socket required for power supply  
**fixed socket 5P 125A 400V IP67**

Cable length up to 50m  
**5 x 35mm<sup>2</sup> (copper cable)  
or 5 x 65mm<sup>2</sup> (aluminum cable)**

External connection protection  
**160A**

### TRAILER DETAILS

<b>Type</b>	<b>Length</b>
<b>LEVEQUES</b>	<b>13.73m</b>
<b>Min. area for placement</b>	<b>Width</b>
<b>15.50m x 6.50m</b>	<b>2.53m</b>
	<b>4.82m (unfolded)</b>
<b>Weight</b>	<b>Height</b>
<b>24 000kg</b>	<b>3.75m</b>

### UNIT FEATURES

- Air-conditioning and recuperation system in all rooms
- Storage cabinets
- Sound system and FM tuner
- Hydraulic patient lift for bed trolley access
- Sink with running water
- IT/Network connections in operator room

### ADDITIONAL FEATURES

- Reference clinician workstation
- CD/DVD burner
- RIS/PACS/Network connection
- Injector: Optivantage Dh - dualhead

### MAINTENANCE AND SERVICE

Maintenance and service contracts included