



GE SIGNA VOYAGER 1.5T MR11

System ID:	MR11
Manufacturer:	General Electric
Model:	Signa Voyager 1.5T
Year of production:	2017
No. of channels:	33
Software:	PX26.0
Gantry opening:	70cm
Max table load:	205kg
Options:	TDI HEAD NECK ARRAY ANTERIOR ARRAY 8ch FOOT ANKLE COIL GEM FLEX INTERFACE 16ch FIXED GE FRU:5430000-23 REV 3 GEM FLEX COIL 16L ARRAY HEAD STAND 16 ch T/R KNEE COIL BREAST COIL HEADREST FLAT TABLE HD FLAT GEM TABLE BREAST ARRAY

ELECTRICITY

Power cable length
15m – 20m

Maximum power parameters
120kW

Type of socket required for power supply
fixed socket 5P 125A 400V IP67

Cable length up to 50m
**5 x 35mm² (copper cable)
or 5 x 65mm² (aluminum cable)**

External connection protection
160A

TRAILER DETAILS

Type	Length
SMIT	13.82m
Min. area for placement	Width
15.50m x 6.50m	2.55m
Weight	Height
31 000kg	4m

UNIT FEATURES

- Air-conditioning and recuperation system in all rooms
- Storage cabinets
- Sound system and FM tuner
- Hydraulic patient lift for bed trolley access
- Sink with running water
- IT/Network connections in operator room

ADDITIONAL FEATURES

- RIS/PACS/Network connection
- Injector

MAINTENANCE AND SERVICE

Maintenance and service contracts included