



GE HDXT 1.5T MR10

System ID:	MR10
Manufacturer:	General Electric
Model:	HDxt 1.5T
Year of production:	2006
Upgrade:	2016
No.of slices:	16
Software:	16.0_03_1638.1
Gantry opening:	60cm
Max table load:	159kg
Coils:	Head, NECK AND SPINE ARRAY 16CH Head Coil CTL Coil Shoulder Coil HD8Ch HD T/R Quad Extremity Coil GP flex Body Array HD 8CH Wrist Array Liberty 9000 8CH Breast Coil HD 8CH NV Array

ELECTRICITY

Power cable length
15m – 20m

Maximum power parameters
120kW

Type of socket required for power supply
Fixed socket 5P 125A 400V IP67

Cable length up to 50m
**5 x 35mm² (copper cable)
or 5 x 65mm² (aluminium cable)**

External connection protection
160A

TRAILER DETAILS

Type	Length
SCHIEN	13.51m
Min. area for placement	Width
15.50m x 6.50m	2.53m
	4.13m (unfolded)
Weight	Height
36 280kg	4m

UNIT FEATURES

- Air-conditioning in all rooms
- Storage cabinets
- Sound system and FM tuner
- Hydraulic patient lift
- Sink with running water
- IT/Network connections in operator room

ADDITIONAL FEATURES

- Reference clinician workstation
- Software: adw 4.4
- CD/DVD burner
- Duplicator
- RIS/PACS/Network connection
- Injector - Spectris Solaris dualhead

MAINTENANCE AND SERVICE

Maintenance and service contracts included